

# KILLS 99.99<sup>%</sup> BACTERIA VIRUSES MOLD

- water-based, non-flammable, non-toxic
- safe for all hard surfaces
- can be used to decontaminate critical/semi-critical medical devices prior to sterilization
- uses unique, water-based nanotechnology

### UNDERSTAND PRODUCT CAUTIONS BEFORE USE.

Always follow all label directions. Store to keep from freezing and out of reach of children. AVOID food contact. AVOID eye/skin contact, inhalation, or ingestion. Do not contaminate food during the use and storage of this product. Rinse surface with potable water before food preparation. This product is **NOT TO BE USED** as a sterilant on any surface/instrument introduced directly into the human body. It may be used to pre-clean/decontaminate critical/ semi-critical medical devices prior to higher level disinfection.

Distributed by Nanotech Sealants Inc. DIN: 02462710 Hamilton, ON.

## DIRECTIONS OF USE:

- 1. Pre-clean surface to remove heavy soil and other contamination.
- 2. Saturate surface until wet, and allow product to penetrate.
- 3. Wait 10 minutes, then use a clean microfiber cloth to wipe product into surface and absorb the excess. Turn cloth frequently.
- 4. When dry, buff surface with clean microfiber cloth to remove any residue and complete sealing of surface.

Use appropriate PPE while applying. Avoid contact with skin/eyes, inhalation, or ingestion.

#### ACTIVE INGREDIENTS:

Alkyl Dimethyl Benzyl Ammonium Chloride ...... 0.05%

Didecyl Dimethyl Ammonium		
Chloride		0.20%

#### FIRST AID:

**Eyes:** Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical attention immediately.

**Skin:** Flush affected area with soap and water. Consult a physician if irritation develops.

**Ingestion:** If this product is swallowed, do not induce vomiting. If vomiting occurs, keep head below hips to prevent breathing liquid into the lungs. Get emergency medical attention immediately.